




EauFlow RO Tanks 2022 Update


As part of our ongoing efforts to improve our customer service, we wanted to make you aware of updates to the technical details for the Eau Flow RO tanks.

We now offer 2 tank sizes:

PRODUCT CODE:	EYST12	EYST22
Water/output capacity L	12	22
Tank dimensions (dia. x h) cm	26 x 44	38 x 58
Tank weight (empty) kg	4	11
Tank weight (full) kg	15	34



FAQ

1) How do I calculate the tank size needed for the practice?	a) No. of dental chairs x bottle size = X b) No. of autoclaves x reservoir size = Y c) X + Y = Z (litres) d) Ensure tank capacity is larger than Z
2) What is maximum recommended tank capacity per RO unit?	34L output (i.e. 1 x 12L and 1 x 22L tanks)
3) How do I join a 12L and 22L tank together?	Using a t-piece (EY20123) 
4) Where can I order a T-piece (EY20123)?	Either through your usual supplier or online: https://cleancert-hygiene.co.uk/product/t-piece-1-4-6-2mm
5) What if the practice needs more water?	We recommend installing a second RO system